

# HABC Gentotype Data Submission Form

Submitting Investigator	
Name	Bernd Meibohm, PhD
Affiliation	University of Tennessee Health Science Center
Phone	(901) 448-1206
Lab Contact	
Lab Contact Phone	
Date Submitted	3/8/05
Polymorphism Information	
Name	Interleukin 8 Receptor, Beta (CXCR2)
RS# or Unique Identifier	rs2230054, rs1126579 and rs1126580
Gene	IL8RB
Chromosome position	2q34-35
Alleles	rs2230054: C785T; rs1126579: T1208C; rs1126580: G1440A
Assay Information	
Genotyping Method	PCR amplification of allele specific primers (PASA)
Genotypes in HWE (Y/N) (attach HWE Form)	Y
Amount of DNA used	10 ng
PCR primers	For rs2230054: FW- TCGTCCTCATCTTCCTGGTC FM- TCGTCCTCATCTTCCTGGTT R- AGTCCATGGCGAAACTTCTG
	For rs1126579: FW- CCATTGTGGTCACAGGATGT FM- CCATTGTGGTCACAGGATGC R- TGCAGAGCTGTCTCACTGGA
	For rs1126580: RW- GTATTTTATAGTAGAGACAGGGTTTGAC RM- GTATTTTATAGTAGAGACAGGGTTTGAT F- CCTCACCCCTTGCCATAAT
PCR Components and Concentrations (Taq, buffer, MgCl <sub>2</sub> , primers, DMSO, other reagents)	Each reaction mixture contained a 1:12,500 dilution SYBR Green I nucleic acid gel stain 10000X in dimethyl sulfoxide (DMSO) (Molecular Probes, Eugene, OR), 0.2 mM dNTP mixture, 200 nM of both forward and reverse primer, 1U Taq DNA Polymerase (Promega, Madison, WI), 6% DMSO, 1X SmartCycler additive reagent (a 5X additive reagent containing BSA at 1 mg/mL, Trehalose at 750 mM, and Tween-20 at 1% v/v) (Cepheid, Sunnyvale, CA), and 10 ng genomic DNA in 1X PCR Buffer (pH 8.3, 10X solution containing 100 mM Tris HCl, 500 mM KCl, and 15 mM MgCl <sub>2</sub> and 0.01% Gelatin) (Sigma, St. Louis, MO).
PCR Cycling conditions (time and temperature for each step, # of cycles, etc.)	The amplification program consisted of initial denaturation of 95°C (5 min) followed by 27 cycles of 95°C (15 s), annealing at 60°C (30 s), and extension at 72°C (45 s). After amplification, melt analysis was performed by heating the reaction mixture from 60°C to 95°C at a rate of 0.2°C/s.
Detection Oligo(s) (if applicable)	For 785 locus, PCR products for sequencing were generated using: Sense primer: 5'-ATGCGGGTCATCTTTGCTGT-3' Antisense primer: 5'-TTGAGGCAGCTGTGAAGGAT-3'
	For 1208 and 1440 locus, PCR products for sequencing were generated using Sense primer: 5'-GGGTTCTCCTTCTCTTCA-3' Antisense primer: 5'-TTACAGGCACTACCACCAC-3'
Other Reaction Conditions (detection reaction components, incubation conditions, gel %, etc.)	The genotyping method was validated using direct sequencing (ABI Prism® 3100, Applied Biosystems, Foster City, CA) after the PCR products were isolated by QIAquick (Qiagen, Valencia, CA). Amplification reactions were also routinely checked for the presence of nonspecific products by 1% agarose gel electrophoresis.
Other Assay Info	Additional info on the assay can be obtained from <i>Clin Chim Acta</i> 341(1-2):93-100.